LuxaPrint Model

DMG Digital Enterprises SE

Elbgaustraße 248 · 22547 Hamburg · Germany · www.dmg-dental.com
092938/2022-04

Instructions for use

English

Product description

Light-curing resin based on (meth)acrylate for the production of dental devices. For DLP / SLA printers that work at wavelengths of 385 nm or 405 nm.

Select material variants according to the wavelength used by the printer.

We recommend using instruments with material parameters that are already defined in the system: e.g. DMG DentaMile Lab5/Lab5 Pro/DMG 3Demax/DMG 3Delite (DMG), D10+/D20+/D30+/D50+ (Rapid Shape), P10+/P20+/P30+/P50+ (Straumann).

Intended purpose

Manufacturing of dental models using 3D printing. The printed model is suitable for producing thermoformed splints.

Limitations of use

LuxaPrint Model is not a medical device according to the EU Medical Device legislation (Directive 93/42/EEC and Regulation (EU) 2017/745).

Intended user

Dentist, dental technician

Notes for use

- During handling (e.g. transferring), the material should only be exposed to daylight / ambient light for as short a time as possible, since lengthy exposure to light can lead to unwanted curing of the material.
- Wear nitrile gloves, protective clothing and goggles when working with the material/printed object.

For optimum properties of the end product:

- Ensure that the instruments used during the production process are in proper condition and comply with machine and material parameters.
- Ensure clean working procedures. Residues on the instruments can cause defects on the printed object.
- Use the material parameters already defined in the system for the cleaning units DMG 3Dewash (DMG), RS wash (Rapid Shape) and P wash (Straumann).
- Allow objects to cure completely.
- Use the material parameters already defined in the system for the light-curing units DMG 3Decure (DMG), RS cure (Rapid Shape) and P cure (Straumann).
- Maximum curing depth during post-curing: For solid objects and when light-curing on both sides, the material can be up to 7 mm thick (at a curing depth of 3.5 mm).
- Thinner layers increase the accuracy of the printed object and increase printing time.
- Position the printed object as flat as possible on the printer's build platform. A flatter positioning increases accuracy and reduces printing time.

Recommended use

Note: Observe the machine and material parameters.

Note: Shake the bottle for at least 1 min before use. However, the LuxaPrint Model Transparent form does not have to be shaken before use.

Note: Maintain a processing temperature of 23 °C ± 2 °C (73.4 °F ± 3.6 °F).

Note: For optimal process reliability, scan the material's RFID code (supported by DMG DentaMile Lab5/Lab5 Pro/DMG 3Demax/DMG 3Delite (DMG), D10+/D20+/D30+/D50+ (Rapid Shape); P10+/P20+/P30+/P50+ (Straumann)).

Finishing the printed object

Production step	
Recommended drip-off time after raising the platform: 10 min	
Note: Finishing should take place immediately after the construction process where possible.	
1.	Carefully remove the printed object from the build platform. When using the printers DMG 3Delite (DMG), D10+ (Rapid Shape) and P10+ (Straumann), start directly at step 2.
2.	Cleaning: Clean the printed object with isopropanol (approx. 99%) using an appropriate cleaning instrument and specified program (e.g. DMG 3Dewash, RS wash or P wash).
2. (Alternatively)	Pre-clean the printed object with ethanol (approx. 96%) or isopropyl alcohol (approx. 99%) for max. 3 min in an ultrasonic cleaner (with a brush if necessary). Optional: Clean openings, holes and gap areas with compressed air. Clean the printed object in a separate container with fresh ethanol (approx. 96%) or isopropyl alcohol (approx. 99%) for max. 2 min.
3.	Post-curing: Post-cure the printed object using the appropriate light-curing unit and the specified program (e.g. DMG 3Decure, RS cure or P cure).
4.	After removing the support, finish the surface of the printed object using the appropriate instruments.

Warnings/precautions

- Keep out of the reach of children!
- Avoid skin contact with the liquid material and components prior to post-curing. In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Avoid breathing in the fumes.
- Follow the manufacturer's instructions for other products that are used with the material(s).
- Observe the safety data sheet!
- Do not use the material where there are known allergies to any of the components or contact allergies.
- Leftover methacrylate-based materials may influence the setting behavior of silicone impression materials.

Product performance characteristics

Cured material	
Flexural strength	≥ 40 [MPa]
Transverse modulus	≥ 1 [GPa]
Shore-D hardness	≥ 50

Uncured material	
Viscosity	≤ 2.5 [Pa*s]
Light sensitivity to fluorescent lamps	≤ 5:00 [min]

Storage/disposal

- Store in a dry place at room temperature (15 – 25°C/59 – 77°F) and protected from light!
- Even low exposure to light can trigger polymerization.
- Always keep the container tightly sealed; carefully seal immediately after every use.
- Do not use after the expiration date!
- The cured material can be discarded together with household waste following consultation with the disposal company and in compliance with regulations.
- Incinerate the resin mixture as hazardous waste in accordance with official regulations.

Composition

Beige: UDMA, TEDMA, CaCO₃, additive.

Grey, Ivory, Transparent: UDMA, TEDMA, TPO, additive. Contains: TPO

Note: A glossary of all symbols is displayed at: <http://www.dmg-america.com/literature>

Caution: US Federal Law restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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