



NEW!



DMG DentaMile Cure RS-XL:

Efficient, high-volume curing



Cure **36 models** in
under **10 minutes**

A smile ahead together



Revolutionize your digital 3D printing workflow with our cutting-edge DentaMile Cure RS-XL, designed specifically for the demands of large-scale dental laboratories.

DMG DentaMile Cure RS-XL

Featuring proven 360° illumination and advanced vacuum technology, the DentaMile Cure RS-XL ensures high-quality curing results every time. Plus, the vacuum technology eliminates the need for gas during curing, significantly reducing operating costs compared to conventional devices.

Features & Benefits

- **360° illumination** for uniform curing in one cycle
- Large chamber for multiple device curing (**up to 36 models**)
- **Vacuum curing** for an oxygen-inhibition-free process
- High-power LEDs for **quick curing cycles (6-10 min.)**
- Intelligent connectivity with DentaMile Lab 5 Pro and 3Demax printers
- Validated workflow for consistent results

The DentaMile Cure RS-XL is the perfect companion to the DentaMile Lab 5 Pro automated 3D printer. Elevate your digital workflow with our bigger, better, and bolder curing unit. **Discover the difference today.**



Size (W/D/H)	25.59 x 31.49 x 41.33 in
Volume	18.50 x 14.76 x 11.02 in.
Curing Time	6-10 min. - depending on material
Connections	WLAN, LAN
Power Supply	230V or 110V
Weight	330.69 lbs



Learn more about the **DentaMile Cure RS-XL**

Ordering Information

Item #	Product Name	Product Contents
171420	DMG DentaMile Cure RS-XL	1 DMG DentaMile Cure RS-XL, 1 starter kit (consisting of 10 optical cleaning cloth, 1 UV protective goggles, 1-5m ethernet cable, 1 USB Stick, corrugated cardboard)

DMG America, LLC
65 Challenger Road, Suite 340, Ridgefield Park, NJ 07660
tel: 800.662.6383 email: info@dmg-america.com
web: www.dmg-america.com

