

LUXAPRINT COMPATIBLE PRINTERS

DMG is dedicated to continuously providing innovative dental materials and solutions that will help drive growth for your practice or lab. Our LuxaPrint suite of resins have been validated in several printers.



LuxaPrint Model



LuxaPrint Ortho



LuxaPrint Ortho Plus
























LuxaPrint Ortho Flex



LuxaPrint Tray

3D Printer
























Curing Unit

3D Printer	Curing Unit	LuxaPrint Model	LuxaPrint Ortho	LuxaPrint Ortho Plus	LuxaPrint Ortho Flex	LuxaPrint Tray
 DentaMile Desk MC-5	 DentaMile Cure MC	✓	✓	✓	✓	✓
 3D Demax	 3D Decure	✓	✓	✓	✓	✓
	 RS Cure	✓	✓	✓	✓	✓
	 P Cure	✓	✓	✓	✓	✓
	 Otoflash	✓ <i>(Grey Only)</i>	✓	✓	✓	✓
 D20+ by RapidShade	 3D Decure	✓	✓	✓	✓	✓
	 RS Cure	✓	✓	✓	✓	✓
	 P Cure	✓	✓	✓	✓	✓
	 Otoflash	✓ <i>(Grey Only)</i>	✓	✓	✓	✓
 P20+ by Straumann	 3D Decure	✓	✓	✓	✓	✓
	 RS Cure	✓	✓	✓	✓	✓
	 P Cure	✓	✓	✓	✓	✓
	 Otoflash	✓ <i>(Grey Only)</i>	✓	✓	✓	✓
 D10+ by RapidShade	 3D Decure	✓	✓	✓	✓	✓
	 RS Cure	✓	✓	✓	✓	✓
	 P Cure	✓	✓	✓	✓	✓
	 Otoflash	✓ <i>(Grey Only)</i>	✓	✓	✓	✓
 P10+ by Straumann	 3D Decure	✓	✓	✓	✓	✓
	 RS Cure	✓	✓	✓	✓	✓
	 P Cure	✓	✓	✓	✓	✓
	 Otoflash	✓ <i>(Grey Only)</i>	✓	✓	✓	✓

LUXAPRINT

COMPATIBLE PRINTERS



Printer Device Name	Curing Unit	LuxaPrint Model	LuxaPrint Ortho	LuxaPrint Ortho Plus	LuxaPrint Ortho Flex	LuxaPrint Tray
 D30+ by RapidShape	 3Decure	✓	✓	✓	✓	✓
	 RS Cure	✓	✓	✓	✓	✓
	 P Cure	✓	✓	✓	✓	✓
	 Otoflash	✓ <i>(Grey Only)</i>	✓	✓	✓	✓
 DentaMile Lab 5 & DentaMile Lab 5 Pro	 3Decure	✓	✓	✓	✓	✓
	 RS Cure	✓	✓	✓	✓	✓
	 P Cure	✓	✓	✓	✓	✓
	 Otoflash	✓ <i>(Grey & Beige Only)</i>	✓	✓	✓	✓
 D50+ by RapidShade	 3Decure	✓	✓	✓	✓	✓
	 RS Cure	✓	✓	✓	✓	✓
	 P Cure	✓	✓	✓	✓	✓
	 Otoflash	✓ <i>(Grey & Beige Only)</i>	✓	✓	✓	✓
 P50+ by Straumann	 3Decure	✓	✓	✓	✓	✓
	 RS Cure	✓	✓	✓	✓	✓
	 P Cure	✓	✓	✓	✓	✓
	 Otoflash	✓ <i>(Grey, Beige & Trans. Only)</i>	✓	✓	✓	✓
 Max Max X Pro 4K by Asiga (385 nm)	 Otoflash <i>(no vacuum or Nitrogen)</i>	✓	✓	✓	✓	✓
 AccuFab-L4D	 FabCure2	✓ <i>(Grey, Ivory, Trans. Only)</i>				
 AccuFab-D1s	 FabCure2	✓ <i>(Grey, Ivory Only)</i>				
 Ackuretta SOL	 Curie	✓ <i>(Grey, Trans. Only)</i>				✓